

EXPERIMENT - 9.

AIM :- Spectrophotometric determination of iron (II) with 1,10-phenanthroline

THEORY :- Iron (II) reacts with 1,10-phenanthroline to form an orange red complex $[(C_6H_5N_2)_3Fe^{+2}]$. The colour intensity is independent of the acidity in pH range 2-9. If iron (III) is present, it can be reduced with hydroxylamine hydrochloride. The absorbance is measured at a wavelength of 515 nm. The variation in the conc. of a given coloured soln. changes the intensity of the transmitted light. The change in light intensity is measured by the instrument called photo colourimeter / colourimeter. When monochromatic light falls on a soln. sample, some light is absorbed and the intensity of the transmitted light is decreased. The decrease in intensity of light is prop. to the thickness of the absorbing medium and the conc. of soln. This may be expressed by Lambert-beers law -

$$A = \log (I_0/I) = \log (I_0/I) = E \cdot b \cdot c$$

where c is the conc. of the soln. expressed in mol/l and 'E' is a constant characteristic of the solute and the wavelength of light, 'E' is called the molar extinction coefficient. 'A' is absorbance or optical density 'b' of solution, b is path length and related to transmittance ($T = I/I_0$)

Yashwini Khanna

Teacher's Signature _____

EXPERIMENT - 9.


AIM :- Spectrophotometric determination of iron(II) with 1, 10-phenanthroline.

APPARATUS :- Burette, volumetric flasks (50ml), cuvettes, funnel, burette stand and colourimeter.

CHEMICALS Mohr's salt solution (ferrous ammonium sulphate;

REQUIRED : $\text{FeSO}_4 (\text{NH}_4)_2 \text{SO}_4 \cdot 6\text{H}_2\text{O}$, 1,10-phenanthroline, hydroxylamine hydrochloride, acetic acid - sodium acetate buffer of pH 4.5 and sulphuric acid (H_2SO_4).

CHEMICAL REACTIONS :-



n = number of phenanthroline molecules reacting with Fe^{2+} .

Jashawini Khanna.

PROCEDURE :-A) PREPARATION OF SAMPLES :-

1. Take 50ml volumetric flask and add 0, 1, 2, 3, 4, 5 ml of FAS soln. in each flask. Let's name them as K, L, M, N, O, P.
2. Then add 2ml of 1, 10-phenanthroline solution to each of these volumetric flasks.
3. Now dilute each volumetric flasks with deionised water to afford a total volume of 50ml. Stop and mix the contents well by shaking for few minutes. Allow the soln. to stand for 10 minutes.
4. The first flask to which 0 ml of FAS is added (ie. no Fe^{+2}) will serve as a blank.

B) TO DETERMINE λ_{MAX}

1. Get two cuvettes issued from lab.
2. Fill one of them with blank soln. (K) and another one of samples containing Fe. Let's say soln. P.
3. Light of single wavelength can be produced by selecting the filter on photocalorimeter usually the ranges goes from 410 nm - 700 nm.
4. Set the filter to 410 nm. Place the cuvette with blank solution, K, in the sample holder.
5. Set the absorbance to 0%.
6. Now place the second cuvette in the sample holder. Measure the absorbance of the solution. Now you have the absorbance at 410 nm for soln. P.

Yashwini Khanna.

Teacher's Signature _____

OBSERVATIONS :-

(i) ABSORBANCE OF THE SOLUTION AT THE HIGHEST CONCENTRATION (10×10^5) AT VARIOUS WAVELENGTHS.

S.NO	WAVELENGTH	ABSORBANCE
1.	400	0.110
2.	420	0.19
3.	470	0.17
4.	500	0.20
5.	530	0.14
6.	620	0.01
7.	660	0.00

(ii) ABSORBANCE OF THE SOLUTION AT DIFFERENT CONCENTRATION AT λ_{MAX} .

S.NO.	CONCENTRATION(N)	ABSORBANCE
K	0	0
L	1.9×10^{-5}	0.04
M	3.9×10^{-5}	0.09
N	5.9×10^{-5}	0.13
O	7.9×10^{-5}	0.17
P	9.9×10^{-5}	0.20
X	unknown	0.12

Yashwini Khanu.

7. By changing the filter to next wavelength each time, repeat steps 4-6. You need to set the absorbance to zero with K everytime you change λ with filter.
8. Now, you have absorbance of soln. P, over a range of wavelength from 410 nm - 700 nm. You will notice that graph b/w absorbance and wavelength takes an inverse parabola shape with λ_{max} around 500 or 480 nm.

c) MEASUREMENT FOR ABSORBANCE FOR SOLUTIONS L to P AT λ_{max} .


1. Set filter paper to λ_{max} obtained at B.
2. Set absorbance to zero.
3. Measure absorbance for solutions L to P now at λ_{max} . Don't disturb filter in between.
4. Now measure absorbance of unknown sample provided to you.
5. Plot absorbance vs concentration for samples L to P. Connecting max points draw a straight line passing through origin.

Yashwini Khanna.

Teacher's Signature _____


(A) TO FIND λ_{MAX} -

ABSORBANCE
VS
WAVELENGTH GRAPH



(B) TO FIND UNKNOWN CONCENTRATION -

ABSORBANCE VS CONCENTRATION



Yashaswini Khanna

6. Using absorbance value for unknown find out its concentration.

GRAPH (OBSERVATION) :-

The conc(X) can be determined from its absorbance value from the graph. A soln. with Fe concentration of $2 \times 10^{-5} \text{ N}$ contains 10 μg of iron. Hence the content in μg can be calculated.

RESULT :- Fe content in the unknown sample is 29.5 μg of Fe.

PRECAUTIONS :-

1. Remove the funnel after fitting the burette with FAS solution.
2. Stopper the flasks and mix the contents well by shaking vigorously for few minutes.
3. Set the absorbance to zero with blank (K) everytime you change the wavelength with filter.
4. Place the cuvette in one particular direction only each time.

Yashaswini Khanna

Teacher's Signature _____

RESULT:-

The Fe content in the unknown sample is

ug of Fe 29.5 ug of Fe.

(ie., $\frac{10}{2 \times 10^{-5}} \times 5.9 \times 10^{-5}$) ug of Fe.

Yashwini Khanna